



The world leader in serving science

ARL ADVANT'X Intellipower™ Series

Analysis of Polymers

Advances in X-ray Fluorescence

ARL ADVANT'X Intellipower Series




ARL ADVANT'X Intellipower Series:

- A new generation of sequential XRF instruments for the 21st century
- Four spectrometers with a choice of power for various levels of performance and speed of analysis:
 - ARL ADVANT'X Intellipower 1200/2500
 - 1200W or 2500W
 - no external water chiller required
 - ARL ADVANT'X Intellipower 3600/4200
 - 3.6 kW – high performance
 - 4.2 kW – top performance


ARL ADVANT'X Family : Features & Benefits

- Built around unique ARL universal goniometer updated for faster operation
- New X-ray tubes (4GN, 5GN) bring higher sensitivity
- Modular approach for flexibility in price / performance
- Multiple sample loading capabilities
- Unique features for highest performance on ultra-light elements analysis


ARL ADVANT'X Intellipower Series: Side View




X-ray tube evolution : closer coupling for higher sensitivities



ARL 8400




ARL 9400



ARL ADVANT'X Series


X-ray tube evolution : thinner windows

- Thinner window for enhanced excitation of ultra-light elements (Be to Cl)



5GN X-ray tube with 50µm Be window (optional)

- 50µm Be window brings enhanced sensitivity compared to 75µm Be window



Exclusive Gearless Goniometer



View of the latest Moiré fringe gearless goniometer, version F45

ARL ADVANT'X Intellipower Series

... 10 Year Warranty
on the angular positioning
of the XRF goniometer

eXtra Performance in sequential XRF

Exclusive Gearless Goniometer

- 5th Generation Goniometer
- Quantitative analysis of any element from Oxygen to Uranium (in basic configuration)
- Elements from Beryllium to Nitrogen when specific multilayer crystals are fitted
- Full capability for analysis of non-routine samples
 - Semi-quantitative analysis (QuantAS)
 - Standardless analysis (UniQuant 5)

Exclusive Gearless Goniometer


- Moiré Fringe positioning
 - Digital precision $< \pm 0.0002^\circ$
- Slewing speed up to 4800 deg/min
- Accuracy ± 0.01 degree on LiF crystals
- Up to 4 collimators
- Up to 9 crystals
- Two detectors with angular range of 0° to 155°
- Fifth generation system with 25 years experience and 2500 gearless goniometers installed worldwide
- Automatic alignment of crystal and detector
 - No human errors & no radiation hazard
- Automatic self quality control

Exclusive Gearless Goniometer

Moiré Fringe Positioning

- Crystal and detector systems are each linked to a circular glass diffraction grating
- Reader heads contain the smaller angled grating and the optical detector
- Crystal grating is rotated by a θ drive, while the detectors are moved at angle 2θ in order to satisfy Bragg's law :
$$n \lambda = 2d \sin \theta$$

- n = order of radiation
- λ = wavelength of radiation
- $2d$ = interplanar spacing of crystal
- θ = angle of incidence on the crystal



Goniometer Principle of Operation

Choice of Collimators, Crystals & Detectors

FEATURE	Sensitivity	Resolution	Elements												
K Spectra			Be	B	C	N	O - F - Na	Mg	Al	Si	P - S - Cl	K	Ca...Ti...Fe	Co...Zn...Sr	
L Spectra			Sn...Yb	Hf...U											
COLLIMAT.															
X-Coarse 2.6°	V.High	Low													
Coarse 0.6°	High	Fair													
Medium 0.25°	Good	Good													
Fine 0.15°	Low	High													
CRYSTAL															
AXBeB	High	Low	Be												
AX20	Good	Low		B											
AX16	High	Low			C										
AX09	Good	Low				N									
AX06	High	Low					O								
AX03	Good	Good					O								
TIAP	Good	Good					O								
ADP	Low	High						Mg							
PET	Good	Fair							Al						
InSb	High	High								Si					
Ge111	Good	High									P				
LiF 200	Good	Good										K			
LiF 220	Fair	High											Ti		
LiF420	Low	V.High													
DETECTOR															
FPC														Zn	
Scintillation															

FPC = flow proportional counter

ARL ADVANT'X Intellipower Series

Comprehensive sample loading systems

Single position



12 positions




Large X-Y auto-sampler




ARL ADVANT'X : Applications


PETROCHEMICALS: Process, Products & Analysis



Sample preparation: Conventional Solid Samples




**Fused Glass Bead
or Polymer disc**
(prepared by hot pressing)




Bulk Sample
(Metal, Glass or Pressed
Powder)

Other Sample Types



Loose Granules

- loss of precision
- loss of accuracy




Liquids


- (oils, solvents)
- no dilution required

Analysis of traces in polymers

- Catalysts for polymerization
 - Al, Ti and Zn
- Charges and additives
 - Plasticizers
 - Lubricants
 - Stabilizing agents
 - Neutralizers
 - Anti-oxidants
 - Pigments
 - Mg, Si, P, S, Cl and Ca



Analysis of traces in polymers



TYPICAL LOD ON ARL ADVANT'X INTELLIPOWER SERIES

ELEMENT	4200W (3 SIGMA) [PPM]	3600W (3 SIGMA) [PPM]	2500W (3 SIGMA) [PPM]	1200W (3 SIGMA) [PPM]
Mg	0.86	0.93	1.11	1.61
Al	0.23	0.25	0.30	0.43
P	0.16	0.17	0.21	0.30
Cl	0.30	0.32	0.39	0.56
Ca	0.14	0.15	0.18	0.26
Ti	0.10	0.11	0.13	0.19
Cr	0.11	0.12	0.14	0.21
Fe	0.07	0.08	0.09	0.13

- Limits of detection for various elements in polymers (100s counting time)


Demands on Ti analysis in polymers

- Ti used as catalyst
- Undesirable rest in the polymer
- Monitoring required
- ± 0.1 ppm on a range from 0.5 to 5 ppm



Repeatability for Ti analysis in polymers

- 10 runs – 60 seconds analysis time
- Sample comes out of spectrometer after each run
- Standard deviation : 0.052 ppm



Repeatability for trace analysis in polymers

RUN	Al [PPM]	Ca [PPM]	Fe [PPM]	Mg [PPM]	Ti [PPM]	P [PPM]	Ce [PPM]
1	71.3	104.4	11.4	60.1	2.1	15.6	24.2
2	72.0	104.2	11.6	59.9	2.0	15.6	24.0
3	71.6	104.1	11.6	60.8	2.0	15.6	23.0
4	71.5	105.1	11.4	60.6	2.1	15.7	23.5
5	72.5	104.9	11.6	61.3	2.0	15.6	23.2
6	73.5	105.2	11.4	61.5	2.0	15.9	23.7
7	72.7	105.6	11.6	61.7	2.2	16.1	23.6
8	73.4	105.8	11.6	60.5	2.1	16.2	22.4
Avg.	72.3	104.9	11.5	60.8	2.0	15.8	23.5
SD	1.4	1.0	0.1	0.2 8	0.1	0.4	1.3

- 8 runs on the same polymer disc
- 20 seconds analysis time per element
- Power: 4200W

References in polymers analysis



Masterbatches

- Schulman Plastics (Belgium)
- Clariant (Italy)

Tapes and adhesives

- MacTac Europe (Belgium)
- Nitto Europe (Belgium)

Polymers films

- Toray Plastics (France)
- Huntsman Chemicals (USA)
- Mitsubishi Chemicals (ex-Hoechst Diafoil) USA

Tyres and rubber

- Pirelli (Italy)

References in polymers analysis



Polyolefins

- Montell (Italy)
- Borealis (Finland)
- Tecnip Qatar (QR)
- SABIC Ibn Zahr (Saudi Arabia)
- Saudi Polyolefins (Saudi Arabia)
- The Polyolefin Company (Singapore)
- Technimont (UAE)
- PolyOne (USA)

Monomers, Polymers, Resins

- Arkema (ex-Atofina Chemicals) (USA)
- Lubrizol (USA)
- BP Amoco Chemicals (USA)
- G.E. Plastics (USA)
- etc..

State of the art in WDXRF: Standard-less analysis


WDXRF has the distinction of being unique and powerful technique by virtue of

Standard-less analysis

Made possible by advancement of mathematical formulation of XRF through Fundamental Parameters and thorough understanding of the physics behind XRF.

Standard-less XRF analysis

Two methods used: Peak hopping vs. Scanning



Versatility of standard-less program UniQuant v.5

Variety of samples types and nature

- Solids such as polymers, alloys and pressed powders
- Liquids with variable mass (catch weight)
- Fused beads with varying dilution (catch weights)
- Small amounts, as little as 10 mg powder.
- Odd shapes, small components.
- Monolayer on substrate. Composition and thickness
- Multilayer on substrate. Thickness of layers.

Limits of detection in oils

Element	Crystal/ Detector	LoD [ppm]	LoD [ppm]	LoD [ppm]
		Intellipower	Intellipower	Intellipower
		3600 W	2500 W	1200 W
Mg	AX06/FPC	2.2	2.65	3.74
Al	PET/FPC	0.67	0.81	1.14
Si	PET/FPC	0.38	0.46	0.65
S	GE111/FPC	0.25	0.30	0.43
Ca	LIF200/FPC	0.17	0.21	0.29
Cr	LIF200/FPC	0.14	0.17	0.24
Mn	LIF200/FPC	0.14	0.17	0.24
Fe	LIF200/FPC	0.15	0.18	0.26
Cu	LIF200/SC	0.12	0.14	0.2
Zn	LIF200/SC	0.11	0.13	0.19
Sn	LIF200/FPC	0.57	0.69	0.97
Pb	LIF200/SC	0.14	0.17	0.24

3 sigma – 100s counting time – 6 micron PP film

ThermoFisher
S C I E N T I F I C

The world leader in serving science

ARL OPTIM'X

High performance at affordable price

Unique WDXRF system

ARL OPTIM'X

- Optimized configuration for optimum analytical results



ARL OPTIM'X: Unique WDXRF product

WD-XRF
SmartGonio™
Flexible, versatile,
elemental analysis




WD-XRF
Multichromators™
Rapid, dedicated
routine analysis

ARL OPTIM'X


■ Side view

- 1 Sample
- 2 X-ray Tube
- 3 High voltage Power Supply
- 4 Low voltage power supply
- 5 Vacuum Pump
- 6 Electronic Rack
- 7 X-ray Tube cooling
- 8 Carte XMI
- 9 SmartGonio™
- 10 Vacuum tank





Close coupling

- Mini X-ray tube (50W - air cooled)
- Close coupling of optics for optimized sensitivity equivalent to 200 W




Superior Spectral Resolution




- 10 times better resolution compared to EDX
- 8 to 10 times higher count rates

ARL OPTIM'X Exclusive SmartGonio™



View of the Moiré fringe SmartGonio™

ARL OPTIM'X Exclusive SmartGonio™




ARL OPTIM'X Exclusive SmartGonio™

Up to 3 crystals fitted !!

Gives access to all elements
from F to U
(when 2 detectors fitted)

Integrated Analytical Assistant



Analytical information on all analyzed elements
Ease of use and help through Analytical Assistant

ARL Optim'X: Unique WDXRF Spectrometer

- Attractive price/performance ratio
- No water cooling required
- Low power consumption
- Very low detector gas consumption
- Most modern and intelligent electronics platform
- Minimized operating costs
 - no external water supply necessary (neither city water nor external water cooler circuit required)
- Modem connection for remote diagnostics

Typical applications

**Oils, gasoline and other
petrochemical products**



Elemental analysis of oils and gasoline

- ISO 14596 or prEN ISO 20884 (S)
- ASTM D2622 (S)
- ISO 15597 (Cl and Br)
- ISO 14597 (Ni and V)
- ASTM D5059 (Pb)
- DIN 51577 (Cl)
- DIN 51391 (Ca, Zn)
- DIN 51431 (Mg)
- DIN 51363 (P)
-



Elemental Analysis of oils and gasoline


SULPHUR IN GASOLINE AND DIESEL

- Environmental regulations worldwide are getting tighter
- Diesel fuel with $S \leq 10\text{ppm}$ is already produced in some European countries
- By 2006 the limit for diesel has been set at 15ppm in USA
- In order to be sure to have less than 10ppm S in the market, the producer must produce at 4ppm S maximum.
- The XRF instrument must have a maximum LoD of 1.5 ppm to do the job

ARL Optim'X: S in oil and gasoline

ASTM D 2622 oder prEN ISO 20884

Limit of detection: 1.4 ppm in 100s
 1.0 ppm in 200s



ARL Optim'X: S in oil and gasoline

Short term repeatability (120 s)

Two samples - new cell at each run

Sample	25ppm	100ppm
Run 1	25.4	100.6
Run 2	25.7	100.7
Run 3	26.4	101.8
Run 4	26	102.3
Run 5	25	98.4
Run 6	25.9	100.5
Run 7	26.7	101.5
Average	25.9	100.8
Std. Dev.	0.58	1.27

ARL Optim'X: Limits of detection

	Fixed channel
Element	LOD [ppm] 120s
Mg	8
Al	3.1
Si	3.2
P	1.5
S	1.2
Ca	1.7
Fe	0.8
Pb	1

Oil matrix - Fixed channels - 120 s counting time

ARL Optim'X: Limits of detection

		Smart Gonio
Element	SmartGonio configuration	LOD [ppm] 120s
Al	PET/FPC	4.2
Si	PET/FPC	4
P	PET/FPC	2
S	PET/FPC	1.7
K	LIF200/FPC	1.4
Ca	LIF200/FPC	1.5
V	LIF200/FPC	1
Cr	LIF200/FPC	1
Mn	LIF200/FPC	1
Fe	LIF200/FPC	1.1
Ni	LIF200/SC	0.6
Cu	LIF200/SC	0.8
Zn	LIF200/SC	0.6
Pb	LIF200/SC	1.7

Oil matrix
SmartGonio™
120 s counting time

FPC : Flow proportional counter
SC : Scintillation counter

Typical applications

- Polymers and masterbatch products



ARL Optim'X: Application in Polymers

■ Limits of detection in Polymers

Elem.	Line	Crystal	Detector	Collima.	Time	Kv/mA	LOD
F	Ka	AX06	FPC	0.29	120	40/1.25	500 ppm
S	Ka	PET	FPC	0.29	120	40/1.25	4.1 ppm
P	Ka	PET	FPC	0.29	120	40/1.25	4.3 ppm
K	Ka	LIF200	FPC	0.29	120	40/1.25	1.5 ppm
Ti	Ka	LIF200	FPC	0.29	120	40/1.25	1.2 ppm
Cr	Ka	LIF200	FPC	0.29	120	40/1.25	1.1 ppm
Br	Ka	LIF200	FPC	0.29	120	40/1.25	0.9 ppm

ARL Optim'X: Application in Polymers

■ Repeatability in Polymers

	Run	F	Si	P	S	K	Ti	Cr	Br
sample 4	1	0.0098	0.0046	0.0023	0.0463	0.0174	0.00014	0.00002	0.0291
sample 4	2	0.0032	0.0045	0.0022	0.0470	0.0173	0.00010	0.00005	0.0290
sample 4	3	0.0211	0.0047	0.0020	0.0463	0.0175	0.00003	-0.00003	0.0289
sample 4	4	0.0239	0.0047	0.0023	0.0467	0.0172	0.00018	0.00000	0.0291
sample 4	5	-0.0147	0.0047	0.0022	0.0471	0.0177	0.00004	-0.00003	0.0290
sample 4	6	0.0343	0.0047	0.0022	0.0477	0.0172	0.00008	-0.00007	0.0290
sample 4	7	0.0305	0.0045	0.0020	0.0469	0.0173	0.00010	-0.00003	0.0290
sample 4	8	-0.0015	0.0046	0.0023	0.0467	0.0172	0.00001	0.00000	0.0291
sample 4	9	0.0239	0.0048	0.0027	0.0468	0.0172	0.00011	0.00008	0.0290
sample 4	10	0.0220	0.0046	0.0028	0.0472	0.0175	0.00011	0.00006	0.0291
sample 4	11	0.0042	0.0046	0.0024	0.0473	0.0174	0.00008	-0.00001	0.0290
	Ave.	0.0143	0.0046	0.0023	0.0469	0.0174	0.00009	0.00000	0.0290
	SD	0.0151	0.00010	0.00023	0.00042	0.00015	0.00005	0.00005	0.00005

All results given in %