

Specification Sheet

IRXross™

Fourier Transform Infrared Spectrophotometer

The IRXross creates a new concept for infrared spectroscopy. It offers the optimal solution for a new era with diverse application requirements.



Hardware

| | |
|-------------------------|--|
| Interferometer | Michelson interferometer (30° incident angle) Equipped with Advanced Dynamic Alignment system Sealed interferometer |
| Optical system | Single-beam optics |
| Beam splitter | Germanium-coated KBr |
| Light source | High-energy ceramic with 3 years guaranteed |
| Detector | DLATGS detector with temperature control MCT detector (Hg–Cd–Te) with liquid nitrogen cooling (option) |
| Wavenumber range | 7,800 to 350 cm ⁻¹ |
| Resolution | 0.25, 0.5, 1, 2, 4, 8, 16 cm ⁻¹ |
| SN ratio | 55,000:1 or higher (with KBr window) 42,000:1 or higher (with KRS-5 window) (4 cm ⁻¹ resolution, 1 min scan, around 2,100 cm ⁻¹ , peak-to-peak) 8,000:1 or higher (with KBr window, MCT measurement) 6,000:1 or higher (with KRS-5 window, MCT measurement) (4 cm ⁻¹ resolution, mirror speed 9 mm/sec, 120 scans, peak-to-peak) 2,000:1 or higher (with KBr window, rapid scan measurement) 1,500:1 or higher (with KRS-5 window, rapid scan measurement) (16 cm ⁻¹ resolution, mirror speed 40 mm/sec, 1 scan, peak-to-peak) |


| | |
|-----------------------------|---|
| Mirror speed | 4-step selection of 2.0, 2.8, 5.0, or 9.0 mm/sec 10, 20, 30, or 40 mm/sec for Rapid Scan (option) |
| Data sampling | He–Ne laser with 12 months guaranteed |
| Gain | Automatic or manual setting (between 1× and 128×) |
| Sample compartment | Equipped with automatic accessory recognition mechanism W200 × D230 × H170 mm Center focus |
| Dimensions | W505 × D578 × H330 mm |
| Weight | 35 kg |
| Installation site | Temperature: 10°C to 35°C* Humidity: 20 to 70% with no condensation *15°C to 30°C for operation 20 to 90%RH with no condensation when using KRS-5 window (moisture-resistant type) at 10°C to 30°C |
| Power supply voltage | AC 100 to 240 V, 50 to 60 Hz |
| Power consumption | 150 VA (for operation), 8 VA (100 V for standby), 25 VA (240 V for standby) |

LabSolutions™ IR Software

| | |
|---------------------------------|---|
| OS | Windows® 10 Pro 64 bit edition |
| Interface | USB 2.0, USB 3.0 |
| Programs | Postrun, Spectrum, Quantitation, Photometric, Time course (option) |
| Measurement | Spectrum measurement, Continuous measurement, Atmospheric correction measurement, Simple measurement mode |
| Hardware monitor | Self-diagnosis function, Status monitor |
| Manipulation functions | Four arithmetic operations, Normalization, Baseline correction, Smoothing, Derivative, X-axis range change, Data correction, Data points thinning/interpolation, Wavenumber-wavelength conversion, X adjust, Time-temperature conversion, Peak pick, Point pick, Film thickness calculation, Data calculation, Degree of coincidence, Deconvolution, Fourier transform, Kubelka Munk conversion, ATR correction, Advanced ATR correction, Kramers Kronig transformation, Atmosphere correction, Degree of coincidence calculation |
| Manipulation functions (option) | Peak split, 3D recalculation, Spectrum extraction from 3D data |
| Analysis support programs | Contaminant Analysis, Pharma Report, Food Additives Identification, Purity Calculation |
| Search functions | Spectrum search (based on similarity), Peak search, Text search, Combination search, Setting of search conditions, Search of user library and commercial library, Creation of user library, Library of approx. 12,000 spectra of organic compounds, polymers, pharmaceutical products, inorganic compounds, food additives, contaminants, etc. |
| Quantitative functions | Multi-point calibration curve method CLS quantitative method PLS quantitative method (option) Photometrics Recalculation function for quantitative and photometric results |
| Printing functions | Report template creation Printing using report templates Easy printing using the ViewPrint function |

| | |
|--------------------|---|
| Validation program | Complies with Chinese, European, US, and Japanese Pharmacopoeias and ASTM (Please confirm the applicable version.) |
| GLP/GMP support | Tree-structured audit trail function Recording of operation logs and data logs (history) Saving by overwriting the same filename is prohibited |
| Security functions | Coordination with LabSolutions security functions User-group based privilege settings |
| Macro functions | Analysis support program <ul style="list-style-type: none"> • 23 programs for supporting analysis related to contaminant analysis, confirmation test, quantitation and film thickness calculation are included as standard. Easy macro function <ul style="list-style-type: none"> • Collective execution of multiple operations by simply arranging operations in the order of the procedure • Execution possible from the PC Macro platform (option) |
| Optional software | Time course, PLS quantitation, Peak split, 3D processing, Macro platform |
| Optional libraries | Contaminant Library Thermal-Damaged Plastics Library UV-Damaged Plastics Library Sadtler Library STJ Library, etc. |
| File formats | Files of JCAMP-DX, ASCII, IRsolution, HYPER-IR can be loaded and saved. |

Dimensions (PC and printer dimensions are examples.)



Units: mm

IRXross and LabSolutions are trademarks of Shimadzu Corporation and affiliated companies in Japan and/or other countries. Windows is a trademark or a registered trademark of Microsoft Corporation in the United States and/or other countries.



Shimadzu Corporation

www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedures.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.
Company names, products/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation, its subsidiaries or its affiliates, whether or not they are used with trademark symbol "TM" or "®".
Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services, whether or not they are used with trademark symbol "TM" or "®".
Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication.