



### The performance difference between MBJ20 and MBJ30

Difference	MBJ20	MBJ30	MBJ30 Performance Introduction
Memory	21 circular records storage	192 circular records storage	/
Swiping	NO	YES	<ol style="list-style-type: none"> <li>1. It will prompt user to swiping card after finishing one test.</li> <li>2. The corresponding data of bilirubin value, user ID number, and test time can be saved after swiping.</li> </ol>
Uploading Data	NO	YES	<ol style="list-style-type: none"> <li>1. The user ID number, bilirubin value and test time can be checked on computer after uploading data.</li> <li>2. Bilirubin value trend graph also can be checked.</li> </ol>
Battery	Non-Rechargeable	Rechargeable	<ol style="list-style-type: none"> <li>1. 2* rechargeable AA batteries. “Ready” light will turn to yellow and flicker during charging.</li> <li>2. Green light will show the full charge.</li> </ol>

Battery charging heat protection	NO	YES	<ol style="list-style-type: none"> <li>1. Red Ready light will warn if batteries are too hot(&gt;48°C) during charge.</li> <li>2. MBJ30 will stop charging automatically.</li> <li>3. The charge will recover after battery temperature is reduced(&lt;44°C)</li> </ol>
----------------------------------	----	-----	---

مشاوره و خرید:

۰۲۱-۷۵۳۱۱

۰۲۱-۲۲۹۱۹۷۲۰

۰۲۱-۲۲۹۱۹۴۷۴