

Brix & Salinity Refractometer



MASTER-BX/S28M

Cat.No.2484


Brix

Sodium chloride

Dual Scale

Range : Brix 0.0 to 33.0% Range : Sodium chloride 0.0 to 28.0g/100g
 Min. scale : 0.2% Min. scale : 0.2g/100g

Size and weight : 3.2×3.4×20.3cm, 105g



Brix scale

Use the Brix scale on the left side for concentrations of beverages, fruit juices, soups, and more.

Sodium chloride scale

The sodium chloride concentration scale is on the right-hand side.

*Concentrations of sodium chloride dissolved in water only. Not detecting the salinity of a complex solution with other ingredients.

ATAGO CO., LTD.
<http://www.atago.net/> overseas@atago.net

Headquarters: The Front Tower Shiba Koen, 23rd Floor
 2-6-3 Shiba-koen, Minato-ku, Tokyo 105-0011, Japan
 TEL : 81-3-3431-1943 FAX : 81-3-3431-1945



HACCP GMP GLP
 ATAGO products comply with HACCP, GMP, and GLP system standards.

- | | | |
|--|--------------------------------|--|
| ATAGO U.S.A., Inc. | TEL : 1-425-637-2107 | customerservice@atago-usa.com |
| ATAGO INDIA Instruments Pvt. Ltd. | TEL : 91-22-28544915, 40713232 | customerservice@atago-india.com |
| ATAGO THAILAND Co., Ltd. | TEL : 662-982-8718-9 | customerservice@atago-thailand.com |
| ATAGO BRASIL Ltda. | TEL : 55 16 3916-6000 | customerservice@atago-brasil.com |
| ATAGO ITALIA s.r.l. | TEL : 39 2 36557267 | customerservice@atago-italia.com |
| ATAGO CHINA Guangzhou Co., Ltd. | TEL : 86 20-38106065 | info@atago-china.com |
| ATAGO RUSSIA Ltd. | TEL : +7(812)339-20-02 | info@atago-russia.com |

* Specifications and appearance are subject to change without notice.